

STELLAR GL

manual

45° of tilt



- Improved clearance for foot propulsion
- Ideal for hand propulsion with accessible rear wheels
- Knees remain low in tilt – wheel under desks, into vans and access other environments with ease
- Transit-ready option
- Short wheelbase provides excellent maneuverability
- Features adjustable seat height, seat depth and back angle
- Anterior and posterior tilt set-ups available



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Multi-adjustable



Adjustable push bar option



Easily adjustable axle position



Simplified actuator mounting style. Different capacity units fitted for optimum performance.



Improved heel clearance

Seat width (standard) 14" – 22"
Seat depth (adjustable) 16" – 22"
Seat height (adjustable) 13" – 18"
Back angle (adjustable) 90° – 122°
Back height 25" or 31"
Dynamic tilt 45°

Overall width
 12" & 16" rear wheels
 20" & 24" rear wheels

seat width + 9.5"
 seat width + 10.5"

Overall length
Weight capacity
Chair weight
Transit option available
Warranty

30" – 37" (without front rigging)
 250 lbs
 42 lbs
 ANSI/RESNA WC-19
 PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. US Patent # 6 447064 B1

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver, BC
 V5X 4N6

F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com